



Anti-BtR1 LDV (aa 938-1037) polyclonal antibody (DPATB-H81235)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-BtR1 LDV Polyclonal Antibody
Target	Btb1
Immunogen	Fusion protein: NRETGLVYVDYETQGSGEVLDRDGDPTHRIFNLIDNFM GEGEGNRNQNDTEVLVILLDVNDNAPELPPPSSELSWTISE NLKQGVRLPHIFAPDRDEP, corresponding to amino acids 938/1037 of BtR1 LDV from Manduca sexta (tobacco hornworm). NRETGLVYVDYETQGSGEVLDRDGDPTHRIFNLIDNFM GEGEG
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Virus
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Btb1 is a currently uncharacterized human ankyrin repeat protein BAA92655. It has been suggested to bind to CUL3
Keywords	BTK Binding Protein; Btb1; BtR1 LDV; cadherin related protein receptor BT R1