



Anti-Decoquinate polyclonal antibody (DPATB-H81369)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Decoquinate Polyclonal Antibody
Target	Decoquinate
Immunogen	Decoquinate conjugated to BSA.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	cELISA
Format	Liquid
Size	50 µl
Buffer	PBS, 0.88% Sodium chloride
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Decoquate is a coccidiostat used in veterinary medicine.

Keywords

6-Decoxy-7-ethoxy-4-oxo-1H-quinoline-3-carboxylic acid ethyl ester
