



# Anti-Dexamethasone polyclonal antibody (DPATB-H81375)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti-Dexamethasone Polyclonal Antibody
<b>Target</b>	Dexamethasone
<b>Immunogen</b>	Chemical / Small Molecule: Dexamethasone conjugated to BSA.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	N/A
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, FC, ICC, IF
<b>Positive Control</b>	Dexamethasone-thyroglobulin (2µg/ml) conjugate can be used for plate coating.
<b>Format</b>	Liquid
<b>Size</b>	1 ml
<b>Buffer</b>	PBS, pH 7.2
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction**

The drug Dexamethasone is a synthetic corticosteroid used to treat inflammation and autoimmune conditions. It is more potent than naturally occurring cortisol and commonly prescribed for various conditions of arthritis, asthma and cancer.

---

**Keywords**

Aeroseb-D; Aeroseb-Dex; anaflogistico; Aphtasolon; Auxiron; Bisu DS; component of Azium; component of Deronil; component of Dronactin; Corson; Corsone; Cortisumman; Decaderm; Decadron Tablets, Elixir; decalix; Dexamethasone

---