



Anti-Diclazuril polyclonal antibody (DPATB-H81377)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Diclazuril Polyclonal Antibody
Target	Diclazuril
Immunogen	Diclazuril (4-N)-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	EIA
Format	Liquid
Size	50 µl
Buffer	20mM PBS, 150mM Sodium chloride, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Diclazuril is a synthetic compound of the triazines group (2,6-Dichloro- α -(4-chlorophenyl)-4-(4,5- dihydro-3,5-dioxo-1,2,4-triazin-2(3H)-yl)- benzeneacetonitrile) a benzeneacetonitrile derivative used for the control of coccidiosis poultry.

Keywords

Diclazuril; Nuoqiu; Vecoxan; Diclazuril premix 0.5%,1%; CLINACOX; DICLAZURIL
