



# Anti-Dicloxacillin polyclonal antibody (DPATB-H81380)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Sheep Anti-Dicloxacillin Polyclonal Antibody
<b>Target</b>	Dicloxacillin
<b>Immunogen</b>	Dicloxacillin conjugated to BSA.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	N/A
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, cELISA
<b>Format</b>	Liquid
<b>Size</b>	50 µl
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	Dicloxacillin is a narrow-spectrum beta-lactam antibiotic of the penicillin class. It is used to treat infections caused by susceptible Gram-positive bacteria. It is active against beta-lactamase-
---------------------	--

producing organisms such as *Staphylococcus aureus*, which wo

---

**Keywords**

) -5-methyl-4-isoxazolyl)carbonyl)amino)-3,3; 4-thia-1-azabicyclo(3.2.0)heptane-2-carboxyl; 4-thia-1-azabicyclo(3.2.0)heptane-2-carboxylicacid,6-(((3-(2,6-dichlorophenyl; 6-(3-(2,6-dichlorophenyl)-5-methyl-4-isoxazolecarboxamido)penicillanicacidqr; 6-dichl

---