



# Anti-Digitoxin polyclonal antibody (DPATB-H81384)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Rabbit Anti-Digitoxin Polyclonal Antibody
Target	Digitoxin
Immunogen	Digoxin-BSA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	RIA
Format	Liquid
Size	250 µl
Buffer	Whole serum
Preservative	15mM Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

## BACKGROUND

Introduction	Digitoxin is a purified cardiac glycoside similar in structure and function to Digoxin. Unlike Digoxin, Digitoxin is eliminated from the body via the liver and not the kidneys. Both drugs are used to treat various heart conditions. Cardiac glycosides bind
Keywords	3beta,5beta)-3-((o-2,6-dideoxy-beta-d-ribohexopyranosyl-(1-4)-o-2,6-dideoxy-be; (3beta,5beta)-hexopyranosyl)oxy]-14-hydroxy; (3beta,5beta)-y)-14-hydroxy;.4)-2,6-dideoxy-beta-D-ribo-hexopyranosyl)oxy]-14-hydroxy-, (3beta,5beta)-;.4)-O-2,6-dideoxy-beta-D-ri