



# Magic™ Anti-E. cloacae Polyclonal antibody (DPATB-H81390)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Chicken Anti-Enterobacter cloacae Polyclonal Antibody
<b>Target</b>	E. cloacae
<b>Immunogen</b>	Whole formaldehyde killed bacteria (Enterobacter cloacae)
<b>Isotype</b>	IgY
<b>Source/Host</b>	Chicken
<b>Species Reactivity</b>	E. cloacae
<b>Purification</b>	IgY fraction
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	0.9% Sodium chloride
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store at +4°C.

## BACKGROUND

**Introduction**

Enterobacter cloacae is a clinically significant Gram-negative, facultatively-anaerobic, rod-shaped bacterium.

---

**Keywords**

Bacteria; Proteobacteria; Gamma Proteobacteria; Enterobacteriales; Enterobacteriaceae; Enterobacter; Enterobacter cloacae; E. c. subsp. cloacae; E. c. subsp. dissolvens; Bacillus cloacae; Bacterium cloacae; Cloaca cloacae; Aerobacter cloacae; Aerobacter c

---