



# Anti-E. coli O63 Polyclonal antibody (DPATB-H81394)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Rabbit Anti-E. coli O63 Polyclonal Antibody
Target	E. coli O63
Immunogen	Tissue / cell preparation E. coli O63 (heat inactivated).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E. coli
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

## BACKGROUND

**Introduction**

Escherichia coli is a gram negative bacillus that belongs to a larger group of Enterobacteriaceae bacteria that inhabit the gastrointestinal tract. Although usually a harmless resident of the gut, some strains have the potential to cause serious problems, es

**Keywords**

E. coli O63; Bacteria; Proteobacteria; Gammaproteobacteria; Enterobacteriales; Enterobacteriaceae; Escherichia; E. coli; Escherichia coli; B