



# Anti-FLJ39060 (aa 1-34) polyclonal antibody (DPATB-H81472)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti-FLJ39060 Polyclonal Antibody
<b>Target</b>	FLJ39060
<b>Immunogen</b>	Recombinant fragment corresponding to amino acids 1-34 of Human FLJ39060 MVDGRTRTIINDIFFTEPTPEMSSLPVRSHSSLS Run BLAST with Run BLAST with
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P
<b>Positive Control</b>	Human colon tissue.
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	pH: 7.20Preservative: 0.02% Sodium azide59% PBS, 40% Glycerol
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

# BACKGROUND

Introduction	Putative uncharacterized protein FLJ39060
Keywords	FLJ39060