



Anti-Galantide (C-terminal) polyclonal antibody (DPATB-H81491)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Galantide Polyclonal Antibody
Target	Galantide
Immunogen	Gly-Trp-Thr-Leu-Asn-Ser-Ala-Gly-Tyr-Leu-Leu-Gly-Pro-Gln-Gln-Phe-Phe-Gly-Leu-Met-NH ₂ (i.e. GWTLNSAGYLLGPQQFFGLM) conjugated to KLH via a carboxyl group GWTLNSAGYLLGPQQFFGLM Run BLAST with Run BLAST with
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Protein G purified
Conjugate	Unconjugated
Applications	EIA, RIA
Cellular Localization	Secreted
Format	Liquid
Size	100 µg
Buffer	50% Glycerol, PBS, pH 7.5
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Galactose Operon Repressor
Keywords	Galantide; L-Methioninamide, glycyl-L-tryptophyl-L-threonyl-L-leucyl-L-asparaginy-L-seryl-L-alanylglycyl-L-tyrosyl-L-leucyl-L-leucylglycyl-L-prolyl-L-glutaminy-L-glutaminy-L-; phenylalanyl-L-phenylalanylglycyl-L-leucyl-; Galanin (1-13)-Substance P (5-1