



Anti- GM130 Polyclonal antibody (DPATB-H81539)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	GM130
Immunogen	Synthetic peptide conjugated to KLH derived from within residues 750 to the C-terminus of Fruit fly (<i>Drosophila melanogaster</i>) GM130. (Immunogen available as DAG-P2053)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC, IFA, WB
Positive Control	Drosophila S2 Cells
Format	Liquid
Size	100 µg
Buffer	1% BSA, PBS. pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Golgi auto-antigen; probably involved in maintaining cis-Golgi structure.
Keywords	130 kDa cis Golgi matrix protein; 130 kDa cis-Golgi matrix protein; Cis golgi matrix protein GM130; GM130; GM130 autoantigen; GOGA2_HUMAN; GOLGA 2; GOLGA2; Golgi autoantigen; Golgi autoantigen golgin subfamily a 2; Golgi matrix protein GM130; Golgin 95; golgin A2; Golgin subfamily a 2; Golgin subfamily A member 2; Golgin-95; MGC20672; SY11 protein;