



# Anti-H. alvei Polyclonal antibody (DPATB-H81586)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Hafnia alvei (conjugated) Polyclonal Antibody
<b>Target</b>	H. alvei
<b>Immunogen</b>	Hafnia alvei total protein extract.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	H. alvei
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, WB
<b>Format</b>	Liquid
<b>Size</b>	50 µl
<b>Buffer</b>	50% Glycerol, ddH2O
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

Hafnia is the genus of the Enterobacteriaceae family whose only species is the Gram-negative, facultatively anaerobic, rod-shaped bacterium *Hafnia alvei*. *H. alvei* is a commensal of the human gastrointestinal tract and not normally pathogenic, but may cause

---

**Keywords**

Hafnia Alve; Bacteria; Proteobacteria; Gamma proteobacteria; Enterobacteriales; Enterobacteriaceae; Hafnia

---