



Anti-Haptoglobin (native) polyclonal antibody (DPATB-H81632)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Haptoglobin Polyclonal Antibody
Target	Haptoglobin
Immunogen	Full length native Haptoglobin purified from bovine plasma
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ELISA
Cellular Localization	Secreted.
Format	Liquid
Size	150 µg
Buffer	pH: 7.50Preservative: 0.02% Sodium azidePBS, 50% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Haptoglobin combines with free plasma hemoglobin, preventing loss of iron through the kidneys and protecting the kidneys from damage by hemoglobin, while making the hemoglobin accessible to degradative enzymes. Uncleaved haptoglobin, also known as zonulin,
Keywords	BP; Haptoglobin alpha chain; Haptoglobin alpha(1S) beta; Haptoglobin alpha(2FS) beta; Haptoglobin beta chain; Haptoglobin, alpha polypeptide; Haptoglobin, beta polypeptide; HP; HP2 ALPHA2; HP2ALPHA2; HPA1S; HPT; HPT_HUMAN; MGC111141; Zonulin;