



Magic™ Anti-Infectious Bronchitis polyclonal antibody (DPATB-H81747)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Infectious Bronchitis (Massachusetts) Polyclonal Antibody
Target	Infectious Bronchitis
Immunogen	Infectious Bronchitis (Massachusetts).
Isotype	IgY
Source/Host	Chicken
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, ELISA, Neut
Format	Liquid
Size	1 ea
Buffer	Whole serum.
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C (add glycerol to a final volume of 50% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Infectious bronchitis is an acute and highly contagious respiratory disease of chickens. The disease is characterized by respiratory signs including gasping, coughing, sneezing, tracheal rales, and nasal discharge. It is caused by a coronavirus. It is an

Keywords

Infectious Bronchitis
