



# Anti-L. biflexa Polyclonal antibody (DPATB-H81761)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti-Leptospira biflexa Polyclonal Antibody
<b>Target</b>	L. biflexa
<b>Immunogen</b>	Tissue / cell preparation (Whole organism).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	L. biflexa
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Conjugation, ICC, IF
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	0.01M PBS, pH 7.2
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

Leptospirosis, a contagious disease affecting both animals and humans and spread by infection with a bacterial pathogen called *Leptospira*, may result in chronic liver and kidney disease and fatality in the dog. *Leptospira* is a flexible, spiral-shaped, Gra

**Keywords**

*Leptospira biflexa*; Bacteria; Spirochaetes; Spirochaetales; Leptospiraceae; *Leptospira*; *L. biflexa*; *L. australis*; *L. canicola*; *L. grippityphosa*; *L. icterohaemorrhagiae*; *L. pomona*; *L. serjoe*