



Anti-L. pneumophila Polyclonal antibody (DPATB-H81771)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Legionella pneumophila Polyclonal Antibody
Target	L. pneumophila
Immunogen	A whole cell preparation of Legionella pneumophila; ATCC#33152.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	L. pneumophila
Conjugate	Unconjugated
Applications	IHC-Fr
Format	Liquid
Size	250 µl
Buffer	PBS, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	Legionellapneumophila is a flagellated gram negative bacterium found primarily in warmwater environments. Legionella pneumophila causes a type of pneumonia calledLegionnaire disease
---------------------	---

and a milder condition called Pontiac fever. Infection is acquired through

Keywords

Proteobacteria; Gamma Proteobacteria; Legionellales; Legionellaceae; Legionella; L. pneumophila; Legionella pneumophila; Legionella pneumophila serogroup 1; Bacteria; γ-proteobacteria; L. longbeachae; Legionella longbeachae
