



Anti-Ochratoxin A (C-terminal) polyclonal antibody [HRP] (DPATB-H81854)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Ochratoxin A Polyclonal Antibody, HRP-Conjugated
Target	Ochratoxin A
Immunogen	The immunogen is prepared by NHS-activated Ochratoxin A via its carboxyl group coupled to KLH
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	HRP
Applications	cELISA
Cellular Localization	Secreted
Format	Liquid
Size	100 µg
Buffer	PBS, 50% Glycerol
Preservative	None
Storage	Shipped at 4°C. Store at -20°C.

BACKGROUND

Introduction

Ochratoxin A is a polyketide-derived dihydroiso-coumarin moiety linked through the 12-carboxy group to phenylalanine. Ochratoxin A is the most important and most commonly occurring of a structurally related group of compounds. It is produced by some speci

Keywords

OA; Ochratoxin A non-ribosomal peptide synthetase; OchratoxinA; OTA;