



Anti-H. pylori Polyclonal antibody (DPATB-H83141)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Helicobacter pylori Polyclonal Antibody
Target	H. pylori
Immunogen	Purified from heat-treated cells of H. pylori strain, CH20429.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	H. pylori
Conjugate	Unconjugated
Applications	IHC-FoFr, IHC-P
Format	Liquid
Size	500 µl
Buffer	20 mM tris-borate, 150 mM sodium chloride, 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	The spiral shaped bacterium <i>Helicobacter pylori</i> is strongly associated with inflammation of the stomach and is also implicated in the development of gastric malignancy. <i>H. pylori</i> is known to
---------------------	--

cause peptic ulcers and chronic gastritis in human. It is associ

Keywords

H. pylori;
