



Anti-HSV type 1 Polyclonal antibody (DPATB-H83158)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-HSV1, prediluted Polyclonal Antibody
Target	HSV type 1
Immunogen	Whole rabbit corneal cells, infected with herpes simplex virus type 1 (strain Mac Intyre), were solubilized with detergent.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HSV
Purification	Protein A purified
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr
Format	Liquid
Buffer	Inert stabilizer
Preservative	0.05% Sodium Azide
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction	Herpes simplex type 1 (HSV-1) belongs to a family that includes HSV-2, Epstein-Barr virus
---------------------	-------------------------------------------------------------------------------------------

(EBV) and Varicella zoster (chicken pox) virus amongst others. HSV-1 and HSV-3 are extremely difficult to distinguish from each other. Members of this family have a?

Keywords

Herpes simplex virus 1; Herpes Simplex Virus; Herpes Simplex Virus Type 1; HSV 1; Human herpesvirus 1; Human herpesvirus type 1;
