



Anti-Theophylline polyclonal antibody (DPATB-H83192)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|--|
| Product Overview | Sheep Anti-Theophylline Polyclonal Antibody |
| Target | Theophylline |
| Immunogen | Theophylline conjugated to Bovine thyroglobulin. |
| Isotype | IgG |
| Source/Host | Sheep |
| Species Reactivity | Bovine |
| Purification | IgG fraction |
| Conjugate | Unconjugated |
| Applications | ELISA, cELISA |
| Format | Liquid |
| Size | 50 µl |
| Buffer | 0.33% Sodium phosphate, 0.88% Sodium chloride |
| Preservative | 0.09% Sodium Azide |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |

BACKGROUND

Introduction

Theophylline is a colorless crystalline alkaloid, derived from tea leaves or made synthetically, used in medicine especially as a bronchial dilator.

Keywords

dimethylxanthine,
