



Anti-SARS-CoV Spike glycoprotein Polyclonal antibody (DPATB-H83218)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-SARS spike glycoprotein Polyclonal Antibody
Target	SARS-CoV Spike glycoprotein
Immunogen	Synthetic peptide: T TFDDVQAPNY TQHTS, corresponding to amino acids 19-35 of SARS spike glycoprotein. T TFDDVQAPNY TQHTS Run BLAST with Run BLAST with
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	SARS-CoV
Purification	Protein G purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	200 µl
Buffer	0.2% Gelatin, PBS
Preservative	0.05% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

A novel coronavirus has been identified as the causative agent of SARS (Severe Acute Respiratory Syndrome). Coronaviruses are a major cause of upper respiratory diseases in humans. The genomes of these viruses are positive stranded RNA approximately 28 to

Keywords

E2,E2 glycoprotein,Human coronavirus spike glycoprotein,Peplomer protein,S,S glycoprotein,Severe acute respiratory syndrome spike glycoprotein,Severe acute respiratory syndrome virus spike glycoprotein,Spike glycoprotein,VGL2,