



Magic™ Anti-Stanniocalcin 1 polyclonal antibody (DPATB-H83229)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Stanniocalcin 1 Polyclonal Antibody
Target	Stanniocalcin 1
Immunogen	Synthetic peptide conjugated to KLH
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Secreted
Format	Liquid
Size	50 µl
Buffer	PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction

Stanniocalcin 1 is a homodimeric phosphoglycoprotein expressed in a wide variety of tissues, including kidney, prostate, thyroid, bone, and ovary. It is considered to have roles in many physiological processes, including bone development, reproduction, wou

Keywords

Stanniocalcin, Stanniocalcin-1, STC, STC1,

GENE INFORMATION

Entrez Gene ID

[6781](#)

UniProt ID

[P52823](#)