



Anti-THC polyclonal antibody (DPATB-H83231)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-THC Polyclonal Antibody
Target	THC
Immunogen	Tetrahydrocannabinol-KLH
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µl
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Tetrahydrocannabinol (THC) is the active chemical in cannabis and is one of the oldest
--------------	--

hallucinogenic drugs known. The cannabinoids belong to a class of chemicals called terpenoids, meaning terpene like. These compounds occur as essential oils within many

Keywords

Tetrahydrocannabinol,THC,
