



Anti-Tetrahydrofolic Acid polyclonal antibody (DPATB-H83237)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Tetrahydrofolic Acid Polyclonal Antibody
Target	Tetrahydrofolic Acid
Immunogen	Folic acid conjugated to BSA
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA
Cellular Localization	Cell membrane
Format	Liquid
Size	100 µl
Buffer	40% Glycerol, Whole serum
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Folic acid and folate (the anion form) are forms of a water-soluble B vitamin. Folate is necessary for the production and maintenance of new cells. This is especially important during periods of rapid cell division and growth such as infancy and pregnancy.
Keywords	Folate, folic acid, THFA, Vitamin B9,