



Anti-Somatostatin polyclonal antibody (DPATB-H83251)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Somatostatin Polyclonal Antibody
Target	Somatostatin
Immunogen	Chemical / Small Molecule conjugated to KLH
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Protein G purified
Conjugate	Unconjugated
Applications	RIA, ELISA
Cellular Localization	Secreted
Format	Liquid
Size	100 µl
Buffer	PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

The cyclic tetradecapeptide somatostatin is widely distributed throughout the body and is an important regulator of endocrine and nervous system function. It exerts its biologic actions by binding to specific high affinity receptors on the cell surface. S

Keywords

Growth hormone release inhibiting factor, SMST, Somatostatin 14, Somatostatin 28, SST,

GENE INFORMATION

Entrez Gene ID

[6750](#)

UniProt ID

[P61278](#)