



Anti-Vitamin B2 polyclonal antibody (DPATB-H83256)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat Anti-Vitamin B2 Polyclonal Antibody
Target	Vitamin B2
Immunogen	Chemical / Small Molecule, Vitamin B2, conjugated to BSA.
Isotype	IgG
Source/Host	Rat
Species Reactivity	N/A
Purification	Preadsorbed IgG fraction
Conjugate	Unconjugated
Applications	ELISA, IHC-Fr
Cellular Localization	Secreted
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Riboflavin (vitamin B2) is manufactured in the body by the intestinal flora and is easily absorbed, although very small quantities are stored, so there is a constant need for this vitamin. It is required by the body to use oxygen and the metabolism of ami
Keywords	Flavin Mononucleotide,FMN,Riboflavin,