



Anti-phenyl Acetic acid polyclonal antibody (DPATB-H83274)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-phenyl Acetic acid Polyclonal Antibody
Target	phenyl Acetic acid
Immunogen	Chemical / Small Molecule; synthetic phenyl acetic acid conjugated to protein carrier.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, ICC, WB, ICC/IF
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Synthetic phenyl Acetic acid is a white crystalline compound with a honey like odour in low concentrations, used in some perfumes and in penicillin G production. In its pure form it is harmful by ingestion, inhalation or through skin contact and is a severe

Keywords

2 phenylacetic acid, Alpha toluic acid, Alpha tolylic acid, Benzene acetic acid, benzeneacetic acid, Benzyl carboxylic acid, Phenyl ethanoic acid, PPA,