



Anti-E. coli O112ac Polyclonal antibody (DPATB-H83332)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-E. coli O112ac Polyclonal Antibody
Target	E. coli O112ac
Immunogen	Tissue / cell preparation of E. coli O112ac (heat inactivated).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E. coli
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction

Escherichia coli is a group of gram-negative bacilli belonging to the family Enterobacteriaceae and is considered as one of the bacteria in the human intestinal normal microbial flora. Their serological types are determined in combination with somatic antigens.

Keywords

Escherichia coli O112ac;