



# Anti-E. coli O114 Polyclonal antibody (DPATB-H83333)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti-E. coli O114 Polyclonal Antibody
<b>Target</b>	E. coli O114
<b>Immunogen</b>	Tissue / cell preparation of E. coli O114 (heat inactivated).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	E. coli
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Agg
<b>Positive Control</b>	Antigenic suspension.
<b>Format</b>	Liquid
<b>Size</b>	2 ml
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C. Do not freeze.

## BACKGROUND

**Introduction** Escherichia coli are Gram negative bacterium that are commonly found in the lower intestine of warm-blooded organisms (endotherms). Their serological types are determined in combination with somatic antigens (O group: O1-O173) and flagella antigens (H type)

---

**Keywords** Escherichia coli O114;

---