



Anti-E. coli Polyclonal antibody (DPATB-H83358)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-E. coli Polyclonal Antibody
Target	E. coli
Immunogen	Mixture of E. coli serotypes.
Isotype	IgG
Source/Host	Goat
Species Reactivity	E. coli
Purification	IgG fraction
Conjugate	Unconjugated
Applications	Conjugation, ICC, IF
Format	Liquid
Size	500 µl
Buffer	0.01M PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Escherichia coli is a gram negative bacillus that belongs to a larger group of Enterobacteriaceae bacteria that inhabit the gastrointestinal tract. Although usually a harmless resident of the gut, some strains have the potential to cause serious problems, es

Keywords

Escherichia coli;
