



Anti-E. coli O111 Polyclonal antibody (DPATB-H83383)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-E. coli O111 Polyclonal Antibody
Target	E. coli O111
Immunogen	Tissue / cell preparation Escherichia coli O serotype 111 (heat inactivated).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	E. coli
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction

Escherichia coli are Gram negative bacterium that are commonly found in the lower intestine of warm blooded organisms (endotherms). Their serological types are determined in combination with somatic antigens (O group: O1-O173) and flagella antigens (H type)

Keywords

Escherichia coli; Escherichia coli O111;
