



Anti-Vasonatrin peptide polyclonal antibody (DPATB-H83427)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Vasonatrin peptide Polyclonal Antibody
Target	Vasonatrin peptide
Immunogen	Synthetic peptide: GLSKGCFGLKLDRI GMSG LG CNS FRY (Human) conjugated to KLH GLSKGCFGLKLDRI GMSG LG CNS FRY Run BLAST with Run BLAST with
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein G purified
Conjugate	Unconjugated
Applications	EIA, RIA
Format	Liquid
Size	100 µg
Buffer	50% Glycerol, PBS, pH 7.5
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

VNP is a chimera of atrial natriuretic peptide (ANP) and C-type natriuretic peptide (CNP). This synthetic peptide possesses the 22-amino acid structure of CNP, which is a cardiovascular selective peptide of endothelial origin and is structurally related t

Keywords

Vasonatrin,
