



Anti-PGCA (aa 1042-1057) polyclonal antibody (DPATB-H83449)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-PGCA Polyclonal Antibody
Target	PGCA
Immunogen	Synthetic peptide: CVRTYWLLGERGCSTRG , corresponding to C terminal amino acids 1042-1057 of Human PGCA. CVRTYWLLGERGCSTRG Run BLAST with Run BLAST with
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Type I membrane protein.
Positive Control	Brain, kidney and spleen.
Format	Liquid
Size	100 µg
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Cyclic GMP (cGMP) is a key messenger in several signal transduction pathways. The intracellular levels of cGMP are maintained by the activity of opposing enzymes: synthesizing gualyl cyclases (GC) and hydrolyzing phosphodiesterases (PDEs). The synthesiz
Keywords	Particulate Guanylate Cyclase A;