



Anti-Hydroxyprogesterone (C-terminal) polyclonal antibody (DPATB-H83453)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Hydroxyprogesterone Polyclonal Antibody
Target	Hydroxyprogesterone
Immunogen	4-PREGNEN-17-OL-3, 20-DIONE 3-CARBOXYMETHYLOXIME conjugated to BSA
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ELISA, RIA
Cellular Localization	Secreted
Format	Liquid
Size	100 µg
Buffer	50% Glycerol, PBS, pH 7.5
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Hydroxyprogesterone is a progestogen structurally related to progesterone that has been used for recurrent miscarriage and various menstrual disorders. Women who have had a spontaneous preterm delivery are at greatly increased risk for preterm delivery in
Keywords	17 alpha hydroxyprogesterone caproate; 17 hydroxy progesterone; Progesterone 17OH; Progesterone8(17OH);