



Anti-Dihydrotestosterone polyclonal antibody (DPATB-H83454)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Dihydrotestosterone Polyclonal Antibody
Target	Dihydrotestosterone
Immunogen	Dihydrotestosterone conjugated to BSA
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Protein G purified
Conjugate	Unconjugated
Applications	RIA, EIA
Cellular Localization	Secreted
Format	Liquid
Size	200 µg
Buffer	50% Glycerol, PBS, pH 7.5
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Dihydrotestosterone is a biologically active metabolite of the hormone testosterone, formed primarily in the prostate gland, testes, hair follicles, and adrenal glands by the enzyme 5 α -reductase by means of reducing the 4,5 double-bond. Dihydrotestosteron

Keywords

5 α DHT; DHT; INN;
