



# Anti-RSV Polyclonal antibody (DPATB-H83468)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Goat Anti-Respiratory Syncytial Virus Polyclonal Antibody
<b>Target</b>	RSV
<b>Immunogen</b>	Human Respiratory Syncytial Virus isolate.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	RSV
<b>Purification</b>	IgG fraction
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P, ELISA, Neut
<b>Cellular Localization</b>	Virion. Host cytoplasm
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	PBS, 10mg/ml BSA
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

Respiratory syncytial virus (RSV) is a major cause of respiratory illness in young children. RSV infection produces a variety of signs and symptoms involving different areas of the respiratory tract, from the nose to the lungs. RSV is a negative sense, en

---

**Keywords**

RS virus,RSV,

---