



# Anti-HAV Polyclonal antibody (DPATB-H83476)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Goat Anti-Hepatitis A Virus Polyclonal Antibody
<b>Target</b>	HAV
<b>Immunogen</b>	Purified virions of strain HM175.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	HAV
<b>Purification</b>	Protein G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, Conjugation
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	0.01M PBS, pH 7.2
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	Hepatitis A Virus (HAV) is a 28nm nonenveloped, spherical, positive stranded RNA virus, classified within the genus hepatovirus of the picornavirus family and is among the smallest and
---------------------	--

structurally simplest of the RNA animal viruses. A single large polypr

---

**Keywords**

HAV;

---