



Anti-Swine influenza virus H1N1 Nonstructural Protein 1 Polyclonal antibody (DPATB-H83486)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Swine H1N1 Nonstructural Protein 1 Polyclonal Antibody
Target	Swine influenza virus H1N1 Nonstructural Protein 1
Immunogen	A synthetic peptide containing the sequence specific to the novel Swine H1N1 Nonstructural Protein 1 (Genbank accession no. ACP41110).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Swine influenza virus
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction

Influenza A virus is a major public health threat, killing more than 30,000 people per year in the USA. In early 2009, a novel Swine H1N1 Neuraminidase virus was identified in specimens obtained from patients in Mexico and the United States. Influenza A v

Keywords

S-OIV NS1, Swine origin Influenza A (H1N1) Nonstructural Protein 1,