



Anti-HCV Core 1a Antigen Polyclonal antibody (DPATB-H83511)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-Hepatitis C Virus Core 1a Polyclonal Antibody
Target	HCV Core 1a Antigen
Immunogen	Recombinant full length protein: Recombinant Core (genotype 1a)
Isotype	IgG
Source/Host	Goat
Species Reactivity	HCV
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	100 µl
Buffer	10mM PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

HCV (Hepatitis C Virus) viral core protein forms the internal viral coat that encapsidates the genomic RNA and is enveloped in a host cell derived lipid membrane. The hepatitis C virus (HCV) core protein represents the first 177 amino acids of the viral p

Keywords

Capsid protein C; HCV 1a nucleocapsid; HCV core 1a; HCV core protein genotype 1a; Hepatitis C Virus core protein; p21 , p19;
