



Anti- IVa2 Polyclonal antibody (DPATB-H83514)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-IVa2 Polyclonal Antibody
Target	IVa2
Immunogen	Fusion protein: GRRPAALQHQQDQPQAHPGQRAARSAPLHRDPDYADEDPA PVERHDPGPGSGRAPTTAVQRKPPQPAKRGDMLDRDAVEQV TELWDRLELLGQTLKSMPTA, corresponding to amino acids 5/104 of Human adenovirus type6IVa2 GRRPAALQHQQDQPQAHPGQRAARSAPLHRDPDYADEDPA PVERHDPGPGSGRAPTTAVQRKPPQPAKRG
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Adenovirus
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Plasma membrane
Format	Liquid
Size	100 µl
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw

cycles.

BACKGROUND

Introduction

IVa2 is a member of the superfamily of ATPases associated with a wide variety of cellular activities, including membrane fusion, proteolysis, and DNA replication.

Keywords

Maturation protein; Protein IVA2;
