



Magic™ Anti-Laryngotracheitis polyclonal antibody (DPATB-H83530)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Laryngotracheitis Polyclonal Antibody
Target	Laryngotracheitis
Immunogen	Tissue / cell preparation: Infectious Laryngotracheitis virus.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Virus
Purification	>95% by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA, Neut
Format	Liquid
Size	200 µl
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C (add glycerol to a final volume of 50% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Laryngotracheitis, also known as LT, is a respiratory disease of poultry caused by a herpesvirus

(Gallid herpesvirus type 1).

Keywords

Gallid herpesvirus type 1; ILTV; infectious laryngotracheitis virus ILTV;
