



Anti-SARS-CoV type Unknown3 Polyclonal antibody (DPATB-H83544)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-SARS Unknown3 Polyclonal Antibody
Target	SARS-CoV type Unknown3
Immunogen	Fusion protein: LTKKNYSELDDEEPMELDYP, corresponding to amino acids 44/63 of SARS Unknown3. LTKKNYSELDDEEPMELDYP Run BLAST with Run BLAST with
Isotype	IgG
Source/Host	Mouse
Species Reactivity	SARS-CoV
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

SARS is not claimed to have been eradicated (unlike smallpox), as it may still be present in its natural host reservoirs (animal populations) and may return to the human population.

Keywords

SARS Sunk3,
