



Anti-Cytomegalovirus Polyclonal antibody (DPATB-H83553)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Cytomegalovirus Polyclonal Antibody
Target	Cytomegalovirus
Immunogen	Cytomegalovirus.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Cytomegalovirus
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ICC, IF, ELISA
Format	Liquid
Size	100 µl
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Cytomegalovirus is a member of the herpes virus group, which includes herpes simplex virus types 1 and 2, varicella zoster virus (which causes chicken pox), and Epstein Barr virus (which
--------------	--

causes infectious mononucleosis). These viruses share a characterist

Keywords

45 kDa immediate-early protein 2; 55 kDa immediate-early protein 1; CMV; Envelope glycoprotein B; gB; Glycoprotein B; HHV 5; HHV5; Human herpesvirus 5; IE1; IE2; UL122; UL122 regulatory protein IE2; UL123; UL123 regulatory protein IE1; UL55; UL55 envelope
