



# Anti-SARS-CoV Matrix Protein Polyclonal antibody (DPATB-H83581)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Rabbit Anti-SARS M Polyclonal Antibody
Target	SARS-CoV Matrix Protein
Immunogen	Synthetic peptide corresponding to the C-terminus.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	SARS-CoV
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA
Cellular Localization	viral envelope
Format	Liquid
Size	100 µg
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	Store at +4°C.

## BACKGROUND

**Introduction** A novel coronavirus has recently been identified as the causative agent of SARS (Severe Acute Respiratory Syndrome). Coronaviruses are a major cause of upper respiratory diseases in humans. The genomes of these viruses are positive-stranded RNA approximat

---

**Keywords** Matrix glycoprotein,

---