



# Magic™ Anti-Tomato spotted wilt virus Core Protein Polyclonal antibody (DPATB-H83590)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Chicken Anti-Tomato Spotted Wilt Virus Core Protein Polyclonal Antibody
Target	Tomato spotted wilt virus Core Protein
Immunogen	Native Tomato Spotted Wilt Virus (TSWV) particles.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Tomato spotted wilt virus
Purification	IgY fraction
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	125 µg
Buffer	0.15M Sodium Chloride, 50mM Tris HCl
Preservative	0.05% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

The Tospoviruses are a genus (Tospovirus) of negative RNA virus found within the family Bunyaviridae. They are the sole group of plant infecting viruses in this family, as all other described members of the Bunyaviridae infect animals. The genus takes its

---

**Keywords**

TSWV,

---