



Anti-PIV type 1 Polyclonal antibody (DPATB-H83604)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-Parainfluenza Virus type 1 (Biotin) Polyclonal Antibody
Target	PIV type 1
Immunogen	Cantell Strain.
Isotype	IgG
Source/Host	Goat
Species Reactivity	PIV
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µl
Buffer	0.01M PBS, pH 7.2. This product contains no stabilizing proteins.
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Human parainfluenza viruses (HPIV) were first discovered in the late 1950s. HPIV is genetically and antigenically divided into types 1 to 4. HPIV 1 to HPIV 3 are major causes of lower respiratory infections in infants, young children, the immunocompromise

Keywords

HPIV 1; hPIV1; Human parainfluenza 1 virus; Human parainfluenza virus type; Parainfluenza; Parainfluenza Type 1;
