



Anti-HCV Nonstructural Protein 5a Polyclonal antibody (DPATB-H83607)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Hepatitis C Virus NS5a Polyclonal Antibody
Target	HCV Nonstructural Protein 5a
Immunogen	E.coli expressed recombinant Hepatitis C Virus NS5a antigen
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	HCV
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ELISA, WB, ICC, IF
Format	Liquid
Size	100 µg
Buffer	PBS pH 7.2, 0.01% sodium azide
Preservative	0.01% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

HCV is classified into the genus Hepacivirus of the family Flaviviridae. Like all the members of the family, HCV is an enveloped, single-stranded, positive-sense RNA virus. Its genome (about 9602 nt) is flanked at both termini by conserved, highly structured

Keywords

HCV; HCV Non-structural protein 5A; Non-structural protein 5A; NS5a; p56;