



# Anti-Shigella Polyclonal antibody (DPATB-H83624)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Rabbit Anti-Shigella Polyclonal Antibody
Target	Shigella
Immunogen	Mixture of <i>S. boydii</i> , <i>S. flexneri</i> and <i>S. dysenteriae</i> .
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Shigella
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ICC/IF
Format	Liquid
Size	100 µl
Buffer	0.01M PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction**

Shigella is a genus of Gram-negative, non-motile, non-spore forming rod-shaped bacteria closely related to Escherichia coli and Salmonella. The causative agent of human shigellosis, Shigella also cause disease in other primates, but not in other mammals.

---

**Keywords**

S. boydii, S. dysenteriae, S. flexneri, S. sonnei, Shigella Boydii, Shigella Dysenteriae, Shigella Flexneri, Shigella Sonnei,

---