



Anti-Streptococcus group A polyclonal antibody (DPATB-H83627)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Streptococcus group A Polyclonal Antibody
Target	Streptococcus group A
Immunogen	Streptococci (Group A) cell preparation.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Streptococcus
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ICC/IF
Format	Liquid
Size	100 µl
Buffer	10mg/ml BSA, 0.01M PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Store in the dark. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

The genus Streptococcus is comprised of a wide variety of both pathogenic and commensal gram positive bacteria which are found to inhabit a wide range of hosts, including humans, horses, pigs and cows. Within the host, streptococci are often found to colo

Keywords

Group A strep, Group A streptococci, Group A streptococcus,
