



Anti-H. influenzae type d Polyclonal antibody (DPATB-H83630)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Haemophilus influenzae d Polyclonal Antibody
Target	H. influenzae type d
Immunogen	Tissue / cell preparation of Haemophilus influenzae d (heat inactivated).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	H. influenzae
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction

Haemophilus influenzae, formerly called Pfeiffer's bacillus or Bacillus influenzae, is a non-motile Gram-negative rod first described in 1894 by Richard Pfeiffer during an influenza pandemic. A member of the Pasteurellaceae family, it is generally aerobic

Keywords

H. influenzae d;
